

Other samples and available tests

The laboratory also accepts non-urine specimens which may be useful adjuncts to urine drugs of abuse testing. Plasma/serum ethanol can be measured alongside urine ethanol (enzymatic assay): blood samples should be collected into fluoride-oxalate (grey-top) tubes. Identification of unknown substances (e.g. in tablets, powders, smoking products) can be carried out using LC-HR-MS with an extensive compound spectral library of >1,000 compounds – NB: contact the laboratory directly prior to sending these specimens to discuss specific requirements. For both plasma/serum ethanol assays and for unknown substance identification additional charges apply.

Sample requirements

Random urine sample (10–20 mL) collected into a plain 30 mL universal container. Please ensure containers are tightly sealed, and are sent in a clear plastic bag separate from the accompanying request form. Please also ensure containers are labelled with the full name and date of birth of the client, and the specimen collection date and time.

Request forms

New request forms are available to download from our website:

www.viapath.co.uk/our-tests/urine-drug-screen

Laboratory Testing

All samples are analysed using a combination of automated immunoassay and liquid chromatography with high-resolution mass spectrometry (LC-HR-MS – see the [Urine Drug Screen](#) section of this leaflet). Creatinine is measured as a marker of sample integrity. Individual drugs/metabolites and drug classes are reported as either 'positive' or 'negative' based on EWTDS 'cut-off' concentrations (www.ewdts.org); note that µg/L is equivalent to ng/mL).

Additional specialist tests by enzymatic assay (urine ethanol) and by LC-HR-MS (see the [Additional Tests](#) section of this leaflet) are available on request, and can be carried out on the same specimen as supplied originally. Please check the Viapath website regularly for service updates – further new tests will be added in response to service requirements.

Contact information

Laboratory address

Toxicology Unit, 3rd Floor, Bessemer Wing, King's College Hospital, Denmark Hill, London, SE5 9RS

Telephone: +44(0)203 299 5881

Fax: +44(0)203 299 5888

Web: www.viapath.co.uk/departments-and-laboratories/drugs-of-abuse-laboratory-at-kings

Email: kch-tr.toxicology@nhs.net

Laboratory hours (telephone enquiries):

Mon-Fri, 09:00-17:30

Senior Clinical Scientists

Simon Handley

Sarah Belsey

Consultant Clinical Scientist, Head of Department

Prof Robert J Flanagan

Service Delivery Manager

Dr Kishor Raja

Author: L Couchman, Authorised: 1st July 2017



Information for Service Users

The Viapath urine drugs of abuse screening service has changed in response to the changes that are occurring in the range of substances encountered. The current process of immunoassay screening followed by amphetamine and/or opiate confirmation has been discontinued. Instead users now receive a single report for specific drugs (now including dihydrocodeine) measured using high resolution-mass spectrometry, with at present immunoassay results for cannabinoids and benzodiazepines, but not barbiturates. Cut-offs (limits below which any results are reported as 'negative') are based on European Workplace Drug Testing Society (EWDTS) guidelines for named analytes, where available. In addition, a range of new tests have been introduced that can be requested separately. These separate tests include barbiturates (by immunoassay when requested), which nowadays are very rarely encountered in clinical practice.

Note: Analyses for medicolegal purposes should be arranged in advance with a specialist forensic laboratory in order to comply with chain-of-custody requirements, etc.

Basic guidance for interpretation

The information below is only a brief introduction to the interpretation of urine drug screening results. For any further information, or for more detailed interpretation of results for specific samples, please contact the laboratory directly.

Urine Drug Screen (UDS)

Creatinine

Creatinine is measured in all urines to check sample integrity. Samples with a creatinine concentration < 0.5 mmol/L will be reported as 'extremely dilute sample; sample integrity failed'. Those with a creatinine concentration 0.5-2.0 mmol/L will be flagged as 'dilute sample'.

Opioids

(cut-offs 300 µg/L, except 6-MAM 10 µg/L)

Morphine - Total morphine is measured (i.e. unconjugated 'free' morphine + **morphine-3-glucuronide** and **morphine-6-glucuronide** as **morphine equivalents**). 'Positive' morphine results can arise from the administration of (i) morphine, (ii) codeine, (iii) diamorphine, (iv) heroin (illicit diamorphine), (v) the antitussive pholcodine, or (vi) from the ingestion of poppy seeds.

Codeine - Total codeine is measured (i.e. unconjugated 'free' codeine + **codeine glucuronide**). 'Positive' codeine results can arise from the administration of (i) codeine (e.g. in over-the-counter analgesic preparations), (ii) heroin, or (iii) from the ingestion of poppy seeds.

Dihydrocodeine (DHC) - Total DHC is measured (i.e. unconjugated 'free' DHC + DHC glucuronide). Positive results can arise only from the administration of DHC.

6-Monoacetylmorphine (6-MAM; also known as 6-acetylmorphine, 6-AM) is a specific metabolite of diamorphine. Concentrations of 6-MAM in urine are low (cut-off 10 µg/L), and it has a short plasma half-life (a few hours), and so a 'negative' 6-MAM result does not exclude prior use of diamorphine. On the other hand, the presence of 6-MAM proves recent administration of diamorphine.

Amphetamine, metamphetamine, MDMA

(cut-offs 200 µg/L)

Amphetamine, metamphetamine ('crystal meth'), and 3,4-methylenedioxymetamphetamine (**MDMA**, 'ecstasy') are analysed separately. Small proportions (< 10 %) of metamphetamine are metabolised to amphetamine following oral administration, hence this compound may also be reported as 'positive' after use of the parent compounds. Amphetamine-like compounds are typically detectable for 2-5 days post-administration.

Cocaine (as benzoylecgonine, cut-off 150 µg/L)

Cocaine is rapidly hydrolysed in blood after insufflation, smoking, or ingestion. The major metabolite found in urine is **benzoylecgonine**, which can be detected for 1-3 days following use of cocaine.

Methadone & EDDP (cut-offs 250 µg/L)

Methadone and its urinary metabolite 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (**EDDP**) are measured separately. This is important because detection of parent drug without metabolite suggests addition of methadone directly to the urine sample.

Benzodiazepines

(immunoassay cut-off 300 µg/L)

Analysis of benzodiazepines is complex since there are a large number of drugs available and many have common metabolic pathways. It is therefore often very difficult to identify exactly which compound(s) have been administered. The immunoassay used detects the major benzodiazepines and their (common) metabolites, e.g. diazepam, nordazepam, oxazepam and temazepam, as well as others such as chlordiazepoxide. The detection window for benzodiazepines is typically up to 7 days, but varies depending on the half-life of the drug, and may be longer in chronic users.

Cannabis (immunoassay cut-off 50 µg/L)

Cannabis is detected in urine primarily as its metabolite 11-nor- Δ -9-tetrahydrocannabinolic acid (11-THC-COOH). A 'positive' cannabis result indicates that the client may have taken cannabis within the last 2-10 days, although this period may be longer (up to 30 days) for chronic users. A 'negative' result indicates that the patient has not taken cannabis within the last 2-3 days, although other factors such as the dose, other drugs taken and urine concentration make this time period hard to specify exactly.

Additional Tests

Barbiturates (immunoassay cut-off 300 µg/L)

Barbiturates are rarely encountered in current drugs of abuse testing, though our immunoassay is available on request to detect use of the most common barbiturates.

Buprenorphine & metabolites (cut-offs 5 µg/L)

Buprenorphine and its urinary metabolites **norbuprenorphine**, **buprenorphine glucuronide** and **norbuprenorphine glucuronide** are measured separately. Following buprenorphine administration (sublingual or intravenous), urinary concentrations of unconjugated buprenorphine are very often 'negative' (< 5 µg/L), and only metabolites can be detected. As with methadone, the detection of parent drug without metabolites suggests addition of buprenorphine directly to the sample.

Mephedrone (cut-off 200 µg/L)

Mephedrone ('meow meow', M-CAT) is a synthetic stimulant derived from cathinone. The UK is reportedly the largest consumer of mephedrone in Europe, yet most immunoassay-based 'screening' methods fail to detect this compound.

Ketamine & norketamine

(cut-off for 'total' 50 µg/L)

Ketamine is used clinically as a sedative and anaesthetic, but may be abused because of its hallucinogenic effects. Ketamine and one of its major metabolites, **norketamine**, are measured separately and reported as a 'total' (ketamine + norketamine) concentration.

Tramadol & metabolites

(cut-off for 'total' 200 µg/L)

Tramadol is an opioid analgesic with significant abuse potential. It is metabolised to **O-** and **N-desmethyltramadol**, both of which are measured separately and are reported together with the parent drug to give a 'total tramadol' concentration in urine samples.

Pregabalin and Gabapentin

(cut-off 500 µg/L)

Pregabalin and **Gabapentin** are used to relieve neuropathic pain, and to help control seizures. Both have significant abuse potential, and are requested individually.

Methylphenidate, Ethylphenidate & Ritalinic Acid

(cut-offs 500 µg/L)

Methylphenidate is a psychostimulant prescribed as first-line treatment in Attention Deficit Hyperactivity Disorder (ADHD), however it also abused for its stimulant properties. **Ethylphenidate** is a metabolite of methylphenidate when co-ingested with ethanol, but is also sold as a 'legal high'. Both are metabolised to **Ritalinic acid**, and all are measured separately.

Novel/new pharmaceutical substances (NPS)

Testing for NPS ('designer drugs', 'legal highs') can be undertaken on request. We have an ever-expanding range of reference compounds available and are continually adding to the list of compounds and metabolites that we have detected in biological samples. Please contact the laboratory before sending specimens to discuss any specific requirements.

